



## ARGUS/Dissolution

The Cutting Edge Application  
Software for Pharmaceutical  
Dissolution Testing

- 256 MB RAM recommended minimum
- 200 MB hard disk space recommended minimum (includes space required during installation)
- Intel® Pentium® III or compatible at 800 MHz or better (single core)
- CD-ROM drive (only for installation)
- High-resolution monitor (1024x768 or higher, 256 colors or better)
- Mouse or other pointing device
- Serial port, PCI slot, ethernet port (the number of ports needed depends on the instruments used)
- Microsoft® Windows® 2000, XP or Vista
- Adobe® Acrobat® Reader (to view manual)